


# Customization

PST's instrument trays can be customized with our unique silicone bar positioning system

PST instrument trays can be customized to adapt to a wide variety of individualized configurations. Our uniquely engineered GridLok System™ allows for easy positioning of the special Plug'nSert™ silicone insert bars. PST silicone products are molded from medical grade silicone. PST trays are manufactured with state-of-the-art durable medical grade plastics. PST trays are engineered with a "pebble-like" interior surface to insure proper penetration of steam and other sterilizing agents under the mat.



The PST GridLok System provides maximum flexibility for custom applications. Plug'nSert silicone bars engineered for use with the GridLok System are designed to be inserted onto the tray surface. They simply plug in with a minimum of effort.



1-800-263-3557

www.sacor.com  
12-300 Steelcase Rd. W.  
Markham ON L3R 2W2

**Customization:**

PST Trays can be customized for individual needs.

**Protection:**

PST Trays protect your instruments from costly damage.

**Quality/Durability:**

PST Trays are manufactured with strict quality standards and are engineered to last.

**Versatility:**


PST Trays can be used with a wide variety of instruments and can be individually customized.

**Cost Savings:**

PST's unique Gridlock System™ allows for the use of one tray for many applications.

**Sterility Assurance:**

The raised mat assures proper penetration of sterilizing agents between the mat and the tray surface.



PST Trays are engineered with a unique "pebble-like" interior surface. This raises the mat to assure proper penetration of the sterilizing agent between the mat and the tray surface. Independent laboratory tests have shown that the space between the mat and the surface is effectively penetrated.



Our latch is designed to last! Compare the quality and manufacturing standards of PST engineered trays and see the difference. Mold flow analysis gives PST products the strength and durability to satisfy your requirements.

### Surgical Instrument Trays By PST

| Part #                                      | Description  | Dimensions W x L x D |                   |
|---|--|----------------------|-------------------|
|   |  | Inches               | Centimeters       |
| <b>MINI TRAYS</b>                           |  |                      |                   |
| 1300A                                       | Mini Tray - Base, Lid, & Mat                       | 1.5 x 2.7 x 0.75     | 3.8 x 6.8 x 1.9   |
| 1310A                                       | Mini Tray - Base (deep), Lid, & Mat                | 1.5 x 2.7 x 1.0      | 3.8 x 6.8 x 2.5   |
| 1320A                                       | Mini Tray - Base (deep), Lid, & Mat                | 1.5 x 2.7 x 1.25     | 3.8 x 6.8 x 3.2   |
| <b>LARGE MICROSURGICAL INSTRUMENT TRAYS</b> |  |                      |                   |
| 1500A                                       | Instrument Tray - Base, Lid, & Mat                 | 10 x 15 x 0.75       | 25.4 x 38.1 x 1.9 |
| 1520A                                       | Instrument Tray - Base, Insert Tray, Lid, & 2 Mats | 10 x 15 x 1.5        | 25.4 x 38.1 x 3.8 |
| 1522A                                       | Instrument Tray - Base (deep), Lid, & Mat          | 10 x 15 x 1.5        | 25.4 x 38.1 x 3.8 |
| <b>PHACO TRAYS</b>                          |  |                      |                   |
| 2341A                                       | Phaco Tray - Base, Lid, & 2 Bars / 1 Slot          | 4.0 x 7.5 x 1.5      | 10.1 x 19 x 3.8   |
| 2343A                                       | Phaco Tray - Base, Lid, & 2 Bars / 3 Slot          | 4.0 x 7.5 x 1.5      | 10.1 x 19 x 3.8   |
| 2361A                                       | Phaco Tray - Base, Lid, & 2 Bars / 1 Slot          | 6.0 x 10 x 1.5       | 15.2 x 25.4 x 3.8 |
| 2363A                                       | Phaco Tray - Base, Lid, & 2 Bars / 3 Slot          | 6.0 x 10 x 1.5       | 15.2 x 25.4 x 3.8 |
| <b>DIAMOND KNIFE TRAYS</b>                  |  |                      |                   |
| 2501A                                       | Knife Tray, 1 Slot - Base, Lid, & 2 Bars           | 2.5 x 6.0 x 0.75     | 6.4 x 15.2 x 1.9  |
| 2502A                                       | Knife Tray, 2 Slot - Base, Lid, & 2 Bars           | 2.5 x 6.0 x 0.75     | 6.4 x 15.2 x 1.9  |
| 2503A                                       | Knife Tray, 3 Slot - Base, Lid, & 2 Bars           | 2.5 x 6.0 x 0.75     | 6.4 x 15.2 x 1.9  |
| 2505A                                       | Knife Tray, 5 Slot - Base, Lid, & 2 Bars           | 4.0 x 6.5 x 0.75     | 10.2 x 16.5 x 1.9 |
| 2571A                                       | Knife Tray, 1 Slot - Base, Lid, & 2 Bars           | 2.5 x 7.5 x 0.75     | 6.4 x 19 x 1.9    |
| 2572A                                       | Knife Tray, 2 Slot - Base, Lid, & 2 Bars           | 2.5 x 7.5 x 0.75     | 6.4 x 19 x 1.9    |
| 2573A                                       | Knife Tray, 3 Slot - Base, Lid, & 2 Bars           | 2.5 x 7.5 x 0.75     | 6.4 x 19 x 1.9    |
| <b>MICROSURGICAL INSTRUMENT TRAYS</b>       |  |                      |                   |
| 2600A                                       | Micro Tray - Base, Lid, & Mat                      | 2.5 x 6.0 x 0.75     | 6.4 x 15.2 x 1.9  |
| 2610A                                       | Micro Tray - Base, Lid, & Mat                      | 2.5 x 6.0 x 1.25     | 6.4 x 15.2 x 3.2  |
| 2700A                                       | Micro Tray - Base, Lid, & Mat                      | 2.5 x 7.5 x 0.75     | 6.4 x 19 x 1.9    |
| 4600A                                       | Micro Tray - Base, Lid, & Mat                      | 4.0 x 6.5 x 0.75     | 10.2 x 16.2 x 1.9 |
| 4700A                                       | Micro Tray - Base, Lid, & Mat                      | 4.0 x 7.5 x 0.75     | 10.1 x 19 x 1.9   |
| 4710A                                       | Micro Tray - Base, Lid, & Mat                      | 4.0 x 7.5 x 1.5      | 10.1 x 19 x 3.8   |
| 6100A                                       | Micro Tray - Base, Lid, & Mat                      | 6.0 x 10 x 0.75      | 15.2 x 25.4 x 1.9 |
| 6120A                                       | Micro Tray - Base, Insert Tray, Lid, & 2 Mats      | 6.0 x 10 x 1.5       | 15.2 x 25.4 x 3.8 |
| 6124A                                       | Micro Tray - Base (deep), Lid, & Mat               | 6.0 x 10 x 1.5       | 15.2 x 25.4 x 3.8 |

PST Corp. warrants its trays and components for one year from the date of purchase against failure due to manufacturing defects. PST will replace defective products at no charge to the customer upon prepaid return of the product to PST.



Plastic Sterilizing  
Tray Corporation


1520 East Street  
Pittsfield, MA 01201



1-800-263-3557

www.sacor.com  
12-300 Steelcase Rd. W.  
Markham ON L3R 2W2

**SETTING THE STANDARD FOR QUALITY**



# PST SCOPE TRAYS

**PST Scope Trays provide you with the protection you need for your expensive endoscopes and cameras.**

PST scope trays are available in a variety of popular sizes to satisfy most of your instrument requirements. The trays are designed to withstand sterilization protocols for steam, ethylene oxide, dry heat, or cold sterilizing solutions. All PST instrument trays are manufactured with durable medical grade polymers.

PST instrument trays are ideal for storing, transporting, and protecting surgical medical instruments. Easy customization can be achieved for a variety of special configurations utilizing PST's unique Plug 'n Sert silicone bar system. This valuable feature allows you to place the silicone bars exactly where you need them for the most secure positioning.

For more information or to order



**1-800-263-3557**

**www.sacor.com**  
**12-300 Steelcase Rd. W.**  
**Markham ON L3R 2W2**

The Sacor logo is a stylized white 'S' on a black background. To the right of the logo is a close-up photograph of a human eye with a bright blue iris. The text for contact information is positioned to the right of the eye.

**Customization:**

PST trays can be customized for individual needs.

**Protection:**

PST Trays protect your instruments from costly damage.

**Quality/Durability:**

PST Trays are manufactured with strict quality standards and are engineered to last.

**Versatility:**



PST Trays can be used with a wide variety of instruments and can be individually customized.

**Cost Savings:**

PST's unique Gridlock System™ allows for the use of one tray for many applications.

**Sterility Assurance:**

The raised mat assures proper penetration of sterilizing agents between the mat and the tray surface.



## Surgical Scope Trays By PST

| Part # | Description                    | Dimensions W x L x D (Inches) |
|--------|--------------------------------|-------------------------------|
| 6500A  | Scope Tray - Base, Lid, & Bars | 2.6 x 12 x 1.5                |
| 6501A  | Scope Tray - Base and Lid      | 2.6 x 12 x 1.5                |
| 6502A  | Scope Tray - Base, Lid, & Mat  | 2.6 x 12 x 1.5                |
| 6510A  | Scope Tray - Base, Lid, & Bars | 2.6 x 16 x 1.5                |
| 6511A  | Scope Tray - Base and Lid      | 2.6 x 16 x 1.5                |
| 6512A  | Scope Tray - Base, Lid, & Mat  | 2.6 x 16 x 1.5                |
| 6520A  | Scope Tray - Base, Lid, & Bars | 3.0 x 22 x 1.5                |
| 6521A  | Scope Tray - Base and Lid      | 3.0 x 22 x 1.5                |
| 6522A  | Scope Tray - Base, Lid, & Mat  | 3.0 x 22 x 1.5                |
| 6550A  | Scope Tray - Base, Lid, & Bars | 6.5 x 18 x 1.5                |
| 6551A  | Scope Tray - Base and Lid      | 6.5 x 18 x 1.5                |
| 6522A  | Scope Tray - Base, Lid, & Mat  | 6.5 x 18 x 1.5                |
| 6560A  | Scope Tray - Base, Lid, & Bars | 6.5 x 18 x 3.0                |
| 6561A  | Scope Tray - Base and Lid      | 6.5 x 18 x 3.0                |
| 6562A  | Scope Tray - Base, Lid, & Mat  | 6.5 x 18 x 3.0                |

PST Corp. warrants its trays and components for one year from the date of purchase against failure due to manufacturing defects. PST will replace defective products at no charge to the customer upon prepaid return of the product.

**For more information or to order CALL:**



Plastic Sterilizing  
Tray Corporation

**1-800-263-3557**

[www.sacor.com](http://www.sacor.com)  
12-300 Steelcase Rd. W.  
Markham ON L3R 2W2